

20th

Anniversary
SINCE 2001



aquapick®

Officially recommended products by
Korean Dental Association

2023 AQUAPICK CATALOGUE


*The Leader
in oral care solutions*





aquapick®

Company Profile



Lobby



Book Cafe



Office Area



Office Area



Product Demonstration Room



Product Display Room



VIP Reception Room



Conference Hall



Office Area



Island Kitchen



Cafeteria



Aquapick Award and Certification Display



Message from CEO

Although improved living standards have enhanced awareness about the care of the teeth and gums, proper care is still not easy to implement in everyday life. This is primarily due to the presence of blind spots where brushing alone does not offer a good solution. Even brushing three times a day still fails to completely prevent small and large problems with the gums and teeth effectively. These blind spots are areas where a brush cannot reach. Pieces of food left around these areas become hotbeds for growing plaque and accelerating the formation of tartar, presenting a threat to the health of the teeth and gums.

Aquapick, which specializes in developing and manufacturing oral care systems, has made efforts to protect and improve your oral health by offering simple and easy solutions to oral care, which requires more than brushing. Thanks to these efforts, Aquapick has won the Presidential Citation which was given by the President of Korea. This is given to the best of best companies that have contributed to social and economic development. Aquapick has taken the lead in the oral care market and promises to continue to afford customers greater satisfaction with better quality and service through intensive investment in R&D and continuous product innovation.

Aquapick is not satisfied with today's achievements in establishing itself as a leading company in the global oral care market and will continue to make every effort to fulfill its social responsibilities by expanding the base of the oral care system and developing various products, enabling more people to take care of their oral health with ease.

Furthermore, Aquapick will pursue diverse activities for the enhancement of customer satisfaction including fulfillment of thorough warranty services. Aquapick has grown with your continued encouragement. We are thankful for your continued faith and trust and we will press forward in our efforts to meet and exceed your expectations with the highest quality products and services.

Thank you.

CEO of Aquapick Co., Ltd.

Kevin Lee





aquapick®

The products that we are selling are not at their highest quality. It is because we are making better products at the moment.

We are not the 'yes man.'

We are the 'no man.'

We always say 'no.'

We will always say "no" until you say "yes."

Fulfilling customer satisfaction, we will pay back our gratitude with greater customer service.





GOOD HABIT


Correct habits for strong teeth and healthy gums




Must-have for
your oral health
Aquapick!



Removes food debris while
massaging gums gently



Prevents peri-implantitis



Ideal for orthodontic patients



Takes care of pregnant
women's teeth (before & after)



Prevention of
kids' tooth decay




Remove bad breath

Aquapick makes teeth and gums healthy and strong,
Aquapick makes your teeth and gums healthy and strong with best oral care products.



HEALTH

Oral health is the beginning of general health.



→ If you clean your teeth carelessly,



Poor oral hygiene accelerates aging process and reduces cognitive ability.

- Journal of American gerontology (Denmark Research Team) -



Greater risks for pneumonia, asthma, and laryngopharyngitis.

- February, 2002, Issue of the British Journal on the Periodontal Membrane-



Greater risk (up to 2.8 times greater than healthy people) for cardiac arrest.

- International Dental Research Association -



Greater risk (up to 7.8 times greater than healthy people) for premature birth.

- International Dental Research Association -



Infection or inflammation in the roots of teeth, may increase the risk of diabetes by lowering insulin levels and cause difficulty in controlling blood sugar levels.


- Taiwanese Dental Health Association-

RESEARCH

Aquapick is an oral health care device with clinically verified results.


Comparison of interdental plaque using Aquapick

Clinical Institution: The Korean Academy of Oral Health
(done by Paired-test)
Researcher: Dr. Bae
(School of dentistry, Seoul National University)



Comparison between Aquapick, dental floss and interdental brush.

Clinical Institution: The Korean Academy of Oral Health
(done by one-way analysis of variance)
Standard deviation 0.67~1.21
Researcher: Dr. Bae
(School of dentistry, Seoul National University)



Aquapick is helpful to prevent gum disease.





aquapick[®]

Aquapick is a daily **must-have**
device for modern people.



| AQ-230



| AQ-120



| AQ-350

NEW

Aquapick Oral Irrigator AQ-350



Officially recommended
by the Korean Dental Association

Product Name: Aquapick Oral Irrigator
Model: AQ-350
Color: Pink / Gray
Components: Main Body, Water Tank, Jet Tip 4EA,
Jet Tip Holder 4EA
Voltage: 100-230V~, 50/60Hz
Power Consumption: 24W
Dimension: 150x90x210mm (WxDxH)
Weight: 710g
Water Tank Capacity: 600ml
Pulsations: 2,400 pulses/minute

Jet Tip Color: ● ● ● ●



**| Pulsating stream of water
at 2,400 pulses/minute**

Effective on removal of food
between tooth & gums



**| Convenient handle
mounting with
ergonomic design**

Effective on removal of food
between tooth & gums



| AUTO START

If you pull up the handle, it
**automatically starts
functioning** by activating
the sensor.
(The world's first development)



| 600ml water tank

Sufficiently usable volume




Aquapick Cordless Oral Irrigator AQ-230



Officially recommended
by the Korean Dental Association

Product Name: Aquapick Cordless Oral Irrigator
Model: AQ-230
Components: Main Body, Water Tank, Charger,
Jet Tip 2EA
Voltage: 100-240V~, 50/60Hz
Power Consumption: 2W
Dimension: 60x85x275mm (WxDxH)
Weight: 280g
Water Tank capacity: 200ml
Pulsations: 1,400 pulses/minute
3 Cleaning Modes: Clean / Soft / Massage

Jet Tip Color: ● ●



- Ergonomic jet tips for optimal cleaning

Optimal oral cleaning effect

0.6mm powerful ultra-fine stream
1,400 pulsating streams per minute



- 3 cleaning modes

Clean/Soft/Massage

Available for
customized management



- Conveniently Portable

Optimal size for
carrying around



NEW

Aquapick Cordless Oral Irrigator AQ-208



Product Name: Aquapick Cordless Oral Irrigator
Model: AQ-208
Color: Pink/White
Components: Main Body, Water Tank, USB Cable,
Pouch, Jet Tip 2EA
Voltage: 5V
Power Consumption: 8W
Dimension: Folded - 30x60x128mm (WxDxH)
Spreaded - 30x60x200mm (WxDxH)
Weight: 200g
Water Tank Capacity: 140ml
Pulsations: 1,600 pulses/minute
3 Cleaning Modes: Clean/Soft/Massage

Jet Tip Color: ● ●



• Hygienic care

The wide inlet design and the complete separation of the water tank make it easy to clean



• Jet tip tray storage bin on the body

You can conveniently and hygienically store the jet tip



• Compact size

Water tank specifically designed to save space
Easy to carry anytime, anywhere



Aquapick Cordless Oral Irrigator AQ-205



Product Name: Aquapick Cordless Oral Irrigator
Model: AQ-205
Color: White / Pink / Mint
Components: Main Body, Water Tank,
Charging Cable, Jet Tip 2EA
Voltage: 3.7V
Power Consumption: 5W
Dimension: 40.8x58x281.5mm (WxDxH)
Weight: 369g
Water Tank capacity: 170ml
Pulsations: 1,700 pulses/minute
3 Cleaning Modes: Clean / Soft / Massage

Jet Tip Color: ● ●



- Easy to Carry

There is a **jet tip hook** inside the water tank at the bottom of the main body, so it can be easily carried.



- Convenient magnetic charging

USB type cable and **magnetic charging method** make it more **convenient** and **safer** to use than general charging cables.



- Large capacity battery

It can be used for up to **30 days** without charging when it is fully charged for **4 hours**.





“Improving Oral Health Globally,”

Aquapick has achieved international **approval for safety and quality**



EU certification



Korean certification
for safe electric
appliances



Japanese
certification



European certification for
meeting the approved
electrowave standard



Recommended by
the Korean academy
of oral health



Korean intellectual
Property office



www.aquapick.com

| **Aquapick Co., Ltd.** |

20F, SK V1 Center, 171, Gasan digital 1-ro, Geumcheon-gu, Seoul, Korea | E-mail. info@aquapick.com

Customer Service: +82-1588-2809